

# ZETTLER® Medical® 800

## Patient handset

### Features:

- Call button
- Reassurance LED
- Locating light for the night
- Service call button, e.g. for cafeteria
- Headphones connector for standard stereo headphones (32 ohms)
- Radio/TV function control (volume, station/channel selection)
- Voice functions (hands-free or discreet mode)
- Special plug-in devices facilitate changes in configuration
- Moisture protection for wipe and immersion disinfection
- Membrane keypad with tactile feedback, buttons are easily located by touching
- Compliance with DIN/VDE 0834, parts 1 + 2, environment exposure class 1 (ambient conditions and EMC)
- Different variants for protection classes A + B
- HTM 2015 compliant
- IP 67 version



### Patient handsets

Patient handsets offer safety and convenience for patients and nursing staff alike, and are essential prerequisites for today's hospital routine. Intelligent call systems enable person-to-person communication between care personnel and patients.

For the ZETTLER® Medical® 800 call system, Tyco has developed a multifunctional handset precisely matched to patient needs.

### Quality Design

The patient's handset features an ergonomic design and is easy to operate via tactile feedback buttons. Room functions, e.g. control of room lighting and bed lighting as well as special functions, such as roller blind control, are easily configurable. A locating light helps find the right buttons during the night. Radio station and TV channel selection can be performed via the handset independently for each bed - even if there is no link to a charging function. The headphones connector permits listening without disturbing other room occupants.

When not needed, the handset can be stored in a flexibly installed mounting bracket. For patients with restricted movement, the handset can also be attached to the bed with a clip.

# ZETTLER® Medical® 800 Patient handset

## Ease of use through voice communication

Of course, the patient's handset comes with all essential care functions. Depending on requirements, a second call type can be programmed in addition to the patient call, e.g. service call. The unit is waterproof and therefore easy to clean and disinfect.

A built-in microphone and loudspeaker permits direct communication and information exchange between patient and personnel. The patient decides whether he or she wishes to use the handset in hands-free mode (stored in the mounting bracket) or to speak discreetly (manual mode). The volume is automatically adjusted.

## Wide variety of options

Depending on functional requirements, the patient's handset is available in several versions. The basic unit has a call button, reassurance light, locating light and waterproof housing - all other functions are optional.

## Flexible plug-in devices

ZETTLER® VarioLine plug-in devices matching the patient's handsets are available. Variants for commercial switch boxes and for hospital installation units can be installed. All versions of the patient's handset - with and without voice function - as well as diagnostic devices and perfusors are supported. The plug-in device automatically detects which version has been connected, and does not require additional programming for handset variant changes.

## Stock Items

Item	Description
Handset	VL Patient 's handset, 1 call button
	VL Patient 's handset, 1 call button (IP67)
	VL Patient 's handset, 1 call button, 1 light button
	VL Patient 's handset, 1 call button, 2 light buttons
	VL Patient 's handset, 1 call button, 2 light buttons (IP67)
	VL Patient 's handset, 1 call button, 2 light buttons, 1 service call button
	VL Patient 's handset, 1 call button, 1 service call button, 2 function buttons
	VL Patient 's handset, 1 call button, 2 light buttons, 1 service call button, 2 function buttons
	VL Patient 's handset, 1 call button, 2 light buttons, 1 service call button, radio/TV control, connector for headphone
	VL Patient 's handset, 1 call button, 2 light buttons, 1 service call button, radio/TV control, connector for headphone, with speech
	VL Patient 's handset, 1 call button, 2 light buttons, 1 service call button, with speech
	VL Patient 's handset, 1 call button, 2 light buttons, connector for headphone, radio/TV, with speech
	VL Patient 's handset, 1 call button, 2 light buttons, with speech
Accessories	VL Retaining bracket with magnet and mounting accessory
	VL Blanket clip
Documentation	Cable extension, 3 m, for VarioLine patient ' handsets
	Configuration adapter patient 's handset
	VarioLine - Installation- and Assembly Instructions VarioLine Patient Handsets, Operating Instructions, English

For further information on how ADT can help you, contact your local ADT Fire and Security office or phone free on **0800 010 999** or visit us at [www.adt.co.uk](http://www.adt.co.uk)

**Head Office:** ADT Fire and Security plc, Security House, The Summit, Hanworth Road, Sunbury-on-Thames, Middlesex TW16 5DB.

ADT reserves the right to modify or withdraw any product or service without notice

**tyco** / Fire & Security

CE

